



Floor Mat

EL-ANC-FM/90x60



Description

The Floor Mat triggers the controller or nurse call system, thereby alerting the caregiver when pressure is placed on it.

It can be placed on the side of the bed or chair to monitor the user at risk for falls, when a resident arises, or in a doorway to alert the caregiver when a resident attempts to depart.

Features

Colour: Grey

Size: 90 x 60 cm

Type: Non-slip

Materials: PVC

Options

Mono or stereo plug



Bedside Fall Mat



Description

The Bedside Fall Mat triggers an alarm when it senses any forces on the mat.

It alerts and assists the caregiver to take care of the user during the day or at night.

It can be placed on the floor next to a bed or chair of a user at risk of falls.

This mat is thicker and softer and can therefore absorb the fall and reduce the degree of injury.



Features

Standalone or can connect to existing nurse call system

Colour: blue/grey

Size: 180 x 75 cm

Materials: PVC

Prices effective 1 July 2016 - exclusive of GST - subject to change



Electrotek

Fall Protection

Chair Mat Sensor

EL-ANC-CM/38x25



Description

Chair Mat Sensor detects the significant change of pressure when the user gets out of the chair or wheelchair, and alerts the caregiver by activating the alarm/s.

The slim design enhances the comfort of the product, and is easy to remove and clean.

It also increases the mobility of moving between different chairs / wheelchairs.

Features

Colour: White

Size: 38 x 25.5 cm

Materials: PVC

Bed Mat Sensor

EL-ANC-BM/76x25



Description

Bed Mat Sensor detects the significant change of pressure when the user gets out of the bed, and alerts the caregiver by activating the alarm/s.

The slim design enhances the comfort of the product, and is easy to remove and clean.

Various sizes increase the flexibility of different usage.

Features

Colour: Grey

Size: 76.5 x 25.5 cm

Materials: PVC

Prices effective 1 July 2016 - exclusive of GST - subject to change



Electrotek

Fall Protection

Bed Mat Cover

EL-ANC-BM/CVR



Description

Bed Mat cover adds extra protection for the Bed Mat Sensor.

It is easy to remove and clean.

Features

Colour: White

Size: To fit Bed Mat Sensor

Materials: Cotton

Mat Controller

EL-ANC-MAT/DM



Description

The Mat Controller gives an audible and visual signal when the alarm activates.

It can be used for bed, chair or floor mats as well as toilet seat sensors or chair restraint.

The alarm delay time is adjustable to suit the caregiver.

Features

Colour: white with green border

Size: 145mm x 90mm x 45mm

Power Source: 3 x AA batteries or AC adaptor

Time of Delay Message: 0, 4 or 8 seconds

Alarm Tone Options: Low, medium or high volume

Low Battery Indicator

Prices effective 1 July 2016 - exclusive of GST - subject to change