



## Bed Wet Sensor

EL-ANC-BW/76X10



### Description

The Bed Wet Sensor gives an audible and visible signal when moisture is detected.

This product is designed to inform the caregiver when incontinence needs to be addressed.

The sensor can be placed under the buttock / body of the individual.

It can only be used with the Wet Controller.

### Features

Standalone

Can be connected to any nurse call system

Colour: blue

Size: 76 x 10 cm

Moisture resistant

## Chair Wet Sensor

EL-ANC-CW/38X25



### Description

The Chair Wet Sensor gives an audible and visible signal when moisture is detected.

This product is designed to inform the caregiver when incontinence needs to be addressed.

The pad can be placed under the buttocks of the individual.

It can only be used with the Wet Controller.

### Features

Standalone

Can be connected to any nurse call system

Colour: blue

Size: 38 x 25 cm

Moisture resistant

Prices effective 1 January 2017 - exclusive of GST - subject to change

### Wet Controller

EL-ANC-WET/CTRL



#### Description

The Wet Controller gives an audible and visual signal when moisture is detected.

It is used with the bed and chair Wet Sensors.

#### Features

Standalone

Can be connected to any nurse call system

Power Source: DC 6V or 3 x AAA batteries

Low Battery Indicator

Colour: white

Size: 134mm x 87mm x 48mm



### Gel Cushion

EL-ANC-GC/40X36



#### Description

The Gel Cushion is a bedsore prevention cushion giving triple distribution of body pressure compared to a rubber cushion.

It is designed to prevent bedsores of users sitting for lengthy periods.

#### Features

Lightweight and easy to carry.

Lycra neoprene fabric with gel which is waterproof, anti-bacterial, breathable and easy to clean.

Unlike foam and sponge cushions this gel cushion is extremely durable.

Colour: Black

Size: 40cm x 36cm x 2cm



## Toilet Seat Sensor

EL-ANC-TS/14X7

EL-ANC-TS/14X4



### Description

The Toilet Seat Sensor detects the significant change of pressure when the user leaves the toilet seat. It then alerts the caregiver by activating the alarm.

Requires the *INDUNA* Mat Controller.

### Features

Standalone

Can be connected to any nurse call System.

Easy to remove and clean

Colour: white

Sizes: 14 x 6.6 cm and 14 x 3.8 cm

Materials: PVC

## Chair Restraint

EL-ANC-CR



### Description

The Chair Restraint Seatbelt can prevent and reduce falls.

When the Hook & Loop strap is opened, it will trigger the alarm monitor or the nurse call system, and produce sound to alert the caregiver.

It is suitable for users with low weight.

Requires the *INDUNA* Mat Controller.

### Features

Standalone

Can be connected to any nurse call System.

Colour: Black with yellow

Size: Adjustable to 1.36m

Materials: PVC

Prices effective 1 January 2017 - exclusive of GST - subject to change



Electrotek

## Personal Comfort Ancillaries

### Mat Controller

EL-ANC-MAT/DM



#### Description

The Mat Controller gives an audible and visual signal when the alarm activates.

It can be used for bed, chair and floor mats as well as toilet seat sensors and chair restraint. By adding a double adaptor a pendant can also be used.

The alarm delay time is adjustable to suit the caregiver.

#### Features

Colour: white with green border

Size: 145mm x 90mm x 45mm

Power Source: 3 x AA batteries or AC adaptor

Time of Delay Message: 0, 4 or 8 seconds

Alarm Tone Options: Low, medium or high volume

Low Battery Indicator

Prices effective 1 January 2017 - exclusive of GST - subject to change